
MEDIGENIC®

Medical Keyboard
Quick Start Guide

Connect the keyboard to a computer with the equipped USB cord or with PS2 adapter (not supplied). Simply plug and play.

The keyboard comes with a flashing indicator  to alert users to clean the keyboard. It has been preset to twelve hour intervals. The timed interval may be changed. The keyboard can also be configured to emit an audible alarm when typed upon as a further alert that cleaning is needed. Follow instructions at <http://www.esterline.com/interfacetechologies/Products/StandardProducts/MedigenicKeyboard/SoftwareUtility.aspx>. Note that these changes can only be done with a USB connection. Once changed, the new settings will function on either USB or PS2.

Prior to cleaning the surface of the keyboard, press  to turn the keyboard off. The light near the key will turn off. There is no need to disconnect the keyboard from the computer.

Clean the keyboard with a hospital-approved germicidal spray or wipe. If using a spray, wipe the surface of the keyboard with a cloth. The presence of liquid on the keyboard surface and the wiping motion will be recorded by the system sensors. This will turn off the flashing indicator  and reset the timed interval.

Your mouse should also be cleaned at this time.

After cleaning, push  to turn on the keyboard. Note that if you forget, the keyboard will automatically turn itself on after one minute.

Keyboard backlighting control:

On: Fn + up arrow  +  Off: Fn + down arrow  + 

The keyboard is designed to be flush-mounted to a tray or stand. Four M4 screws (not supplied) for a standard 75mm VESA mount can be used.

There are also adhesive-backed mounting feet supplied with the keyboard. Adhere four small mounting feet to the bottom of the keyboard for normal desktop usage. Two additional sets of rear mounting feet in different sizes are supplied to adjust the keyboard angle for user preference. Simply select the desired set and install over the small rear feet.

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